

imFCS: Correlator

New Correlation Job:

Image File: D:\dkfz\FCSTOOLS\trunk\QuickFit3\plugins\tcspcimporter\assets\testdata\alexa2.t3r

Details: output: file format: PicoQuant TITR File prefix: ./results/ postfix: _corr%correlator%_%COUNTER%

input range: from start 0,000 s to end 10,000 s duration: 300s countrates: [2.0599; 0.588542; 0; 0] kHz à 30 seconds

☒ FCS: cut input into: 5 runs

<input checked="" type="checkbox"/> ACF 1	<input checked="" type="checkbox"/> CCF 1,2	<input type="checkbox"/> CCF 1,3	<input type="checkbox"/> CCF 1,4
<input type="checkbox"/> CCF 2,1	<input checked="" type="checkbox"/> ACF 2	<input type="checkbox"/> CCF 2,3	<input type="checkbox"/> CCF 2,4
<input type="checkbox"/> CCF 3,1	<input type="checkbox"/> CCF 3,2	<input type="checkbox"/> ACF 3	<input type="checkbox"/> CCF 3,4
<input type="checkbox"/> CCF 4,1	<input type="checkbox"/> CCF 4,2	<input type="checkbox"/> CCF 4,3	<input type="checkbox"/> ACF 4

the correlator will calculate $CCF(i,j)$, where channel i will be undelayed and channel j will be delayed.

correlator: bin and correlate: Multi-Tau 1 (one monitor per lag)

Correlators (S): 20 # Lags/Cor. (P): 16 # lag factor (m): 2 min. lag time τ_{min} : 1,00 μ s spanned correlator lags: $\tau_{min} = 1\mu s \dots \tau_{max} = 8.38861s$

☒ countrate: 1000,00 μ s

<input checked="" type="checkbox"/> CH 1	<input checked="" type="checkbox"/> CH 2	<input type="checkbox"/> CH 3	<input type="checkbox"/> CH 4
--	--	-------------------------------	-------------------------------

☒ add to project

+ add job

Progress: Overall Progress: 0% max. parallel jobs: 2

running jobs will be put here ... they are processed in the background, so more jobs may be added. click on „Close“ only when all are finished (green check mark)!

Help Close

input range set to „full“ (both checked)

select a t3r file

general TCSPC file settings

FCS/FCCS evaluation settings

countrate evaluation settings

job management and progress

create 5 „runs“, i.e. cut the input stream into 5 segments and calculate one correlation function for each

calculate the ACF for channels 1 & 2 and calculate the cross-correlation function between these channels

configure correlator

calculate binned countrate with 1ms binning for channels 1 and 2